ClinCal[®] - Calibrator

Serum-Calibrator lyophilised, for Vitamin C

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The values were assigned with the ClinRep[®] complete kit for vitamin C. The uncertainty of the mean value is $\pm 3.3\%$ at a confidence level of 95%. The mean value is traceable to NIST-SRM970.

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Intended use:

ClinCal[®] Serum-Calibrators are used for calibration of the HPLC-system. This lyophilised calibrator has to be prepared like a patient sample.

Reconstitution:

Add exactly 0.5 ml HPLC-water to the vial and mix for max. 3 min. After that the lyophilisate is reconstituted and ready to use.

Caution: Immediately after reconstitution the serum-calibrator has to be prepared according to the instruction manual of the ClinRep[®] complete kit for vitamin C and analysed.

Storage and stability:

Originally closed and stored below -18 °C in the dark, the lyophilised serum-calibrator is stable for 24 months, but not longer than the expiration date printed on the label. After reconstitution, the serum-calibrator has to be used immediately and cannot be stored. The average residual moisture of this lot is 0.84 %.

Pack size:

ClinCal[®] Serum-Calibrator for Vitamin C: 5 x 0.5 ml, **order no.: 28013**

Caution:

The human serum which was used for manufacturing the calibrator was tested for the following infectious markers and found negative: HIV 1/2- and HCV-antibodies, hepatitis B-surface antigen. Nevertheless, the serum-calibrator should be considered as potentially infectious and treated with appropriate care.

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