ClinCal® - Calibrator

Serum-Calibrator lyophilised, for 25-OH Vitamins D2/D3

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The values were assigned with the ClinRep® complete kit for 25-OH Vitamins D2/D3. The uncertainty of the mean values is $\pm 2.1\%$ for 25-OH-Vitamin D3 and $\pm 1.9\%$ for 25-OH-Vitamin D2 at a confidence level of 95%. The mean values are traceable to NIST-SRM972.

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Intended use:

ClinCal® Serum-Calibrators are used for calibration of the HPLC-system. This lyophilised calibrator has to be prepared like a patient sample.

Reconstitution:

Add exactly 2.0 ml HPLC-water to the vial and mix for 15 min. When all material is dissolved, the solution is ready to use.

Storage and stability:

Originally closed and stored below -18 °C in the dark, the lyophilised serum-calibrator is stable for 36 months, but not longer than the expiration date printed on the label. After reconstitution the stability of the analytes is:

- at least 3 days when stored at $15 - 30\text{ °C}$ in the dark
- at least 7 days when stored at $2 - 8\text{ °C}$ in the dark
- at least 3 months when stored below -18 °C in the dark
(avoid repeated freezing and thawing)

The average residual moisture of this lot is 1.98 %.

Pack size:

ClinCal® Serum-Calibrator for 25-OH Vitamins D2/D3: 5 x 2 ml, **order no.: 35013**

Caution:

This ClinCal® Calibrator was manufactured from pooled calfserum. It should be considered as potential infectious and treated with appropriate care.

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