ClinCal® - Calibrator

Plasma Calibrator Set lyophilised / Plasma-Kalibrator-Set lyophilisiert

FOR ANTIMYCOTICS / FÜR ANTIMYKOTIKA

Intended use:

ClinCal® plasma calibrators are used for calibration of the analysis-system. These lyophilised calibrators are based on human plasma and are available as a calibrator set with 5 levels, including a blank. The calibrators have to be prepared like patient samples in one series of analyses.

Reconstitution:

Add exactly 1.0 ml of HPLC-water to the vial and mix for 15 min. When all material is dissolved, the solution is ready to use.

Storage and stability:

Originally closed and stored at 2 - 8 °C, the lyophilised plasma calibrators are stable for 48 months, but not longer than the expiration date printed on the label.

After reconstitution the stability of the analytes is:

- at least 8 hours when stored at 15 - 30 °C
- at least 7 days when stored at 2 - 8 °C
- at least 30 days when stored below -18 °C (avoid repeated freezing and thawing)

Notes:

The mean values of the analytes are chosen in a way that a linear calibration curve can be generated over a suitable range of concentration. According to quality assurance all ClinCal® calibrators have to pass strict quality control procedures during manufacturing. RECIPE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume is max. +/- 0.25 %. The average residual moisture of this lot is 1.02 %

Mean values:

The mean values were established in additional independent reference laboratories with supervision of RECIPE, according to the guideline of the German Medical Council (RiliBäk) with statistic methods.

Pack size:

ClinCal® Plasma Calibrator Set for Antimycotics
Level 0 - 4
5 x 1 x 1 ml, **order no.: 9939**

Precautions:

The human plasma which was used for manufacturing the calibrators was tested for the following infectious markers and found negative: HIV 1/2- and HCV-antibodies, hepatitis B-surface antigen, HIV- and HCV-genome (PCR). Nevertheless, the plasma calibrators should be considered as potentially infectious and treated with appropriate care.

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Level 0

This blank-level is manufactured from the same pool of plasma like the others, but it was not spiked with any analyte. Therefore the concentration of each of the analytes is zero.

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RECIPE CHEMICALS + INSTRUMENTS GmbH