## ClinCal® - Calibrator

### Urine Calibrator Set lyophilised / Urin-Kalibrator-Set lyophilisiert

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##### Intended use:

ClinCal® urine calibrators are used for calibration of your analyzing system. These lyophilised calibrators are based on human urine and are available as a calibrator set with 7 levels, including a blank. The calibrators have to be prepared like patient samples in one series of analyses.

##### Reconstitution:

Add exactly 3.0 ml of HPLC-water to the vial and mix for 15 min. When all material is dissolved, the solution is ready to use.

##### Storage and stability:

Originally closed and stored at 2 - 8 °C, the lyophilised urine calibrators are stable for 36 months, but not longer than the expiration date printed on the label.

After reconstitution the stability of the analytes is:

- at least 12 hours when stored at 15 - 30 °C
- at least 7 days when stored at 2 - 8 °C
- at least 3 months when stored below -18 °C (avoid repeated freezing and thawing)

##### Notes:

The mean values of the analytes are chosen in a way that a linear calibration curve can be generated over a wide range of concentration. According to quality assurance all ClinCal® calibrators have to pass strict quality control procedures during manufacturing. RECIPE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume is max. +/- 0.25 %. The average residual moisture of this lot is 0.51 %.

Creatinine concentrations of the different calibrator levels are between 0.5 g/l and 1.50 g/l.

##### Mean values:

The mean values were established in additional independent reference laboratories with supervision of RECIPE, according to the guideline of the German Medical Council (Rilibäk) with statistic methods.

##### Pack size:

ClinCal® Urine Calibrator Set for Nicotine and -Metabolites Level 0 - 6  
7 x 1 x 3 ml, **order no.: MS3013**

##### Precautions:

These calibrators were manufactured from human urine. Therefore the material has to be considered as potentially infectious and should be treated with appropriate care.

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RECIPE CHEMICALS + INSTRUMENTS GmbH