ClinChek® - Control Plasma Control lyophilised / Kontrollplasma lyophilisiert

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Intended use:

ClinChek® plasma controls are used for internal quality assurance in clinical-chemical laboratories. These lyophilised controls are based on human plasma and are available with mean values in two different ranges of concentration. The controls have to be prepared like patient samples in one series of analyses.

Reconstitution:

Add exactly 0.5 ml of HPLC-water to the vial and mix for 15 min. When all material is dissolved, the solution is ready to use.

If it is necessary to have a defined buffer-constitution of the control for further sample preparation (e.g. for pH-adjustment), you can also dissolve the lyophilised plasma control directly in the appropriate buffer.

Storage and stability:

Originally closed and stored below -18 °C, the lyophilised plasma control is stable for 24 months, but not longer than the expiration date printed on the label.

After reconstitution the stability of the analyte is:

- at least 8 hours when stored at 15 - 30 °C
- at least 4 days when stored at 2 - 8 °C
- at least 21 days when stored below -18 °C (avoid repeated freezing and thawing)

Notes:

Routine analyses of the controls are performed according to the test procedure for patient samples. The concentrations of the analyte are chosen in ranges where valid results can be obtained. According to quality assurance all ClinChek® controls have to pass strict quality control procedures during manufacturing. RECIPE guarantees the same stability and constitution for each vial of one lot. The variation of the filling volume is max. +/- 0.5 %. The average residual moisture of this lot is 1.37 %.

Mean values:

The mean values and confidence intervals were established in additional independent reference laboratories with supervision of RECIPE, according to the guideline of the German Medical Council (RiliBäk) with statistic methods.

Pack size:

ClinChek® Plasma Control

Level I

10 x 0.5 ml, **order no.: 32080**

Level II

10 x 0.5 ml, **order no.: 32081**

Level I, II

2 x 5 x 0.5 ml, **order no.: 32082**

Precautions:

The human plasma which was used for manufacturing the controls was tested for the following infectious markers and found negative: HIV 1/2- and HCV-antibodies, hepatitis B-surface antigen, HIV- and HCV-genome (PCR). Nevertheless, the plasma controls should be considered as potentially infectious and treated with appropriate care.

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RECIPE CHEMICALS + INSTRUMENTS GmbH